

ONE STEP.
ONE MINUTE.
ONE CHOICE.
FOR DAILY DISINFECTION



Why Choose Sani-Cloth® AF³?

Sani-Cloth® AF³ wipes are alcohol free, ready-to-use low linting, disposable 2 in 1 cleaning and disinfectant wipes for non-invasive medical devices; non-porous hard surfaces. Sani-Cloth® AF³ can be used to pre-clean or decontaminate medical devices prior to sterilisation.



	Test Method	Organism	Contact Time
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (d)
	EN 16615	Staphylococcus aureus	60 sec (d)
	EN 16615	Enterococcus hirae	60 sec (d)
	EN 13727	Pseudomonas aeruginosa	60 sec (d)
	EN 13727	Staphylococcus aureus	60 sec (d)
	EN 13727	Enterococcus hirae	60 sec (d)
	EN 13727	MRSA- methicillin-resistant staphylococcus aureus	60 sec (d)
Yeast	EN 16615	Candida albicans	60 sec (d)
	EN 13624	Candida albicans	60 sec (d)
	EN 13624	Candida auris	60 sec (d)
Enveloped virus	EN 14476	Vaccinia Virus	60 sec (d)
	EN 14476	Corona,HIV, Hep B & C, SARS, Flu	60 sec (d)
Non-enveloped virus	EN 14476	Norovirus	60 sec (d)
TB	EN 14348	Mycobacterium terrae	60 sec (d)

Dirty Conditions (d)



Equipment Protection: The largest alcohol-free compatibility database in the world, tested with more than 450 medical devices.



Efficient: 1 minute total contact time including TB & norovirus in dirty conditions.



Skin Integrity: Dermatologically tested and PHMB-free. Can be used without gloves.

Sani-Cloth® AF³ is compatible with our Point of Use dispensers for improved compliance!

Recommended Guidance:

1 pack of EAF3011 or 2 packs of EAF3003



Product Description	Product Code	Wipe Size	Classification	PAO	Shelf life
200ct Soft Pack	EAF3011	267mm x 210mm	MD Class IIa	6 months	24 months
100ct Soft Pack	EAF3003	200mm x 200mm	MD Class IIa	3 months	24 months