

SPONGE COUNTER BAGS



E6500N Sponge counter bags

designed to facilitate and improve the post-procedural counting of surgical sponges.

- Assists and organizes the counting process in the OR.
- Easy open pockets for the simple disposal of used sponges
- Accommodates five Laparotomy sponges or ten 4" (10cm) x 4" (10cm)
- Can be hung from sponge counter rack which attaches to most I.V. poles.
- The clear polyethylene pockets allow sponges to be viewed for visual estimate of patient blood loss
- Additional bags can be placed over filled ones as necessary
- Accurate measurement for assessing absorbed fluid
- Packaged in convenient dispenser box
- Excellent visibility for counting sponges
- Reduces body fluid contact
- Latex free

Optional:
E6558N Sponge Counter Dispenser Rack